

## Session 4.1 – Tipificación basada en perfil alélico o gene-by-gene

Isabel Cuesta

BU-ISCIII

Unidades Comunes Científico Técnicas – SGSAFI-ISCIII

04-15 Noviembre 2019, 2ª Edición  
Programa Formación Continua, ISCIII

# Index

- Typing resolution
- Concepts: homology, core, accessory and pan-genome
- Schema definition
- e.g. *Listeria monocytogenes*

# Typing methods: DNA-based methods

## **PFGE** *Gold standard*


Is a rather time-consuming and labour-intensive technique.

Discriminatory power of PFGE profiles is limited as only nucleotide changes in the restriction enzyme recognition sites are detected.

Relatedness of strains may be over- or underestimated

Epidemiologically unrelated isolates may be assigned to one 'pseudo'-cluster whereas even highly related strains fall into distinct clusters.

# Sequence data for taxonomy and typing



Different levels of sequence information can be associated with different taxonomic levels.

The need for higher-resolution characterization of isolates has led to the development of a wide range of strain-typing methods

# Variability between bacterial genomes of the same species

## GENOME EVOLUTION

**Vertical transfer:** is the passing of genetic material by descent

**Horizontal transfer:** is the movement of genetic material among bacteria that do not necessarily share a mother cell.

- transformation: the uptake of DNA by a cell
- conjugation: transfer facilitated by conjugative elements
- phage-mediated transduction

- Point mutations: SNPs, single nucleotide insertion / deletion.
- Large insertions / deletions
- Genome rearrangements
- Transfer of exogenous DNA
- Plasmids or phages

# Concepts

**Homology:** share a common ancestor, either by descent or recombination. No such thing as “significantly homologous”. A sequence either is or is not homologous.

Infer homology from knowledge of evolutionary relationships and from degrees of similarity between sequences, features or other data.

**Orthologues:** sequences have common ancestor and have split due to speciation event.

**Paralogues:** genes arise by gene duplication


# Concepts

**Core genome:** the number of shared features in a pool of genomes. Shared genes among multiple strains are mostly related to house-keeping genes or central metabolic processes, most of the structural information and main genotypic features. **Orthologues** in all genomes of bacteria belonging to the same taxa

**Accessory genome or adaptative genome:** includes genes conferring adaptive advantages to the strain in order to survive in a specific environment. In most cases, these factors are linked to antibiotic resistance, virulence, capsular serotype, adaptation, and might reflect the organisms predominant lifestyle.

**Pangenome:** The term “pan-genome” refers to pan (from Greek παν, whole) and genome (genome) referring to the inclusion of the core and the dispensable genome.

# General analytical process for cgMLST / wgMLST





# Gene-by-gene: Defining a schema

## Sequence-clustering algorithms:

BDBH:

OMCL


COGtriangles

| name | option  |  |
|------|---------|--|
| BDBH | default | Starting from a reference genome, keep adding genomes stepwise while storing the sequence clusters that result of merging the latest bidirectional best hits, as illustrated in Figure 3.                          |
| COGS | -G      | Merges triangles of inter-genomic symmetrical best matches, as described in PubMed= <a href="#">20439257</a> . Note that a single sequence might occasionally be included in several COGS clusters with option -x. |
| OMCL | -M      | OrthoMCL v1.4, uses the Markov Cluster Algorithm to group sequences, with inflation (-F) controlling cluster granularity, as described in PubMed= <a href="#">12952885</a> .                                       |


## BLAST

- % of coverage in the pairwise alignments query/subject
- % of sequence identity in query/subject pairs
- Genome uses as reference genome

# BDBH




# Get\_Homologues






# Get\_Homologues

Contreras-Moreira, Appl Environ Microbiol 2013



50 Streptococcus proteomes from 14 species  
BLAST: minimum pairwise alignment coverage of 75%



# Reasons why schemas are different

Van Tonder et al., PlosCompBiol 2014

Any collection of isolates is a subset of the entire population for the species of interest, and if the subset of isolates has limited genetic diversity then the number of “core” genes shared by all isolates in that sample will be higher than in a dataset which is genetically more diverse.

More generally, the **size of the core genome is dependent on the size of the data set**, with the core genome decreasing in size as more genomes are added to the analysis


## Technical reason:

**Incomplete or “draft” genomes.** This is acceptable for most studies, but analyses of these genomes may exclude a gene from a list of core genes simply because it contains a sequence gap or is otherwise incomplete at that locus in the assembly of one or a few genomes

## BLAST parameters

# Discriminatory power of typing methods

## *Listeria monocytogenes*



The Simpson's index is used to quantify the probability that two unrelated strains are assigned to different typing groups

# Core-genome schemas for *Listeria monocytogenes*

**Ruppitsch et al., 2015**

SeqSphere+



**Pightling et al, 2015**

Bioinformatics pipeline that takes raw sequence reads as input and calculating a core genome profile by comparing it to an expandable database to compile a phylogeny

**Moura et al., 2016**

1748 loci

All 4  
lineages



## BLASTN

Minimum nucleotide identity 70%  
Alignment length coverage 70%

62% of its CDS


Moura et al., Nature Microbiology 2016



# Moura core-genome schema for *Listeria monocytogenes*


## Validation of the cgMLST scheme with a set of 650 genomes

A




Impact of the N50 assembly size in the number of missing loci. Cut-off N50 of 20kb

B




Distribution of the number of missing loci per genome

A




B




1370 Distinct  
allelic profiles in  
1696 genomes


C




D




E




# Comparison with other genome-based MLST schemes




# General analytical process for cgMLST / wgMLST



# Reproducibility of cgMLST allelic calls, sequencing depth & assembly strategie



# BIGSdb-1m



<http://bigsdb.pasteur.fr/listeria>

# Retrospective validation of whole genome sequencing-enhanced surveillance of listeriosis in Europe, 2010 to 2015

Van Walle et al., eurosurveillance 2018

**2,726 Lm isolates** from human cases from 27 EU/EEA countries, 2010-2015

1,069 isolates -> public health laboratories

1,657 isolates -> commercial sequencing provider

MiSeq (2x150, 2x250, 2x300) NextSeq (2x150), HiSeq (2x100), Ion Torrent PGM

## Trimming

- Removal of any adaptor sequences
- Removal of leading bases and trailing bases with  $Q < 25$
- Window of 20 bases has average  $Q < 25$
- Removal reads with length  $< 36b$

**Assembly:** Spades 3.7.1 or Velvet 1.1.04. Minimum contig length 300nt. Assembly with the highest N50 was retained.

**Allelic profile:** two subsets of isolate pairs,  $AD \leq 7$  (closely related isolates likely to share a common epidemiological link),  $AD \leq 150$  (sublineages where isolates are still likely to have common phenotypic properties that may be relevant, e.g. source attribution)

Moura (Bionumerics: reads (kmers) or BLASTN), 1748 loci

Ruppitsch (SeqSphere +3.4.1), 1701 loci

# Retrospective validation of whole genome sequencing-enhanced surveillance of listeriosis in Europe, 2010 to 2015

Van Walle et al., eurosurveillance 2018

## Conclusions:

- The average coverage up to around 55x before trimming and 45x after trimming (Illumina)
- Assembly-based allele calling outperforms reads-based allele calling -> more loci were detected (increase typeability) and the average distances between isolates were slightly smaller.
- Velvet including k-mer optimisation performed slightly worse than SPAdes, but both produced near-equivalent results.
- cgMLST analysis to share assembled genomes rather than sequence read data (fastq , for SNPs analysis or to verify the analyses in some cases, e.g. multi-country outbreaks)
- Important: individual differences in Ads between Moura and Ruppitsch CG schemes can be relatively large, since only 1,261 loci are common to both schemes.
- The  $AD \leq 7$  cutoff is useful for cluster detection for both schemas, in general, although there are exceptions when there are epidemiologically linked isolates with more allele differences than the cutoff (more than one strain or specific sublineages may have higher average mutation rates -> sublineages specific cutoff ??)
- WGS-enhanced surveillance of listeriosis: many clusters found involved more than one country. Earlier detection of clusters.
- The molecular typing results must also be combined with epidemiological and food exposure investigations.



# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



**Table 1. Number and percentage of laboratories submitting results for each method**

|                            | Serotyping        |                |      |       | Cluster analysis |          |      |       |
|----------------------------|-------------------|----------------|------|-------|------------------|----------|------|-------|
|                            | Conventional only | Molecular only | Both | Total | PFGE-only        | WGS-only | Both | Total |
| Number of participants     | 1                 | 12             | 5    | 18    | 3                | 8        | 4    | 15    |
| Percentage of participants | 6%                | 67%            | 28%  | 90%*  | 20%              | 53%      | 27%  | 75%*  |

*Thirteen of the 20 participants (65%) completed both parts (serotyping and cluster analysis) of the EQA.*

*\* Percentage of total number of participating laboratories (20)*

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing

## Annex 7. Reported sequencing details

| Sequencing performed | Protocol (library prep) | Commercial kit  | Sequencing platform |
|----------------------|-------------------------|---|---------------------|
| In own laboratory    | Commercial kits         | Nextera XT DNA library Preparation Kit*   | HiSeq2500           |
| In own laboratory    | Commercial kits         | NEBNext® Fast DNA Fragmentation & Library Prep Set for Ion Torrent, New England Biolabs** | Ion Torrent PGM     |
| Externally           | Commercial kits         | Illumina  | HiSeq 2500          |
| In own laboratory    | Commercial kits         | Ion Xpress™ Plus Fragment Library Kit for AB Library Builder™ System                      | IonTorrent S5XL     |
| In own laboratory    | Commercial kits         | Nextera XT  | MiSeq               |
| In own laboratory    | Commercial kits         | NEXTERA   | MiSeq               |
| In own laboratory    | Commercial kits         | SureSelect QXT Library Prep Kit (Agilent)   | MiSeq               |
| In own laboratory    | Commercial kits         | Nextera XT DNA Library Preparation Kit  | MiSeq               |
| In own laboratory    | Commercial kits         | Nextera XT  | MiSeq               |
| In own laboratory    | Commercial kits         | Nextera XT***   | Miniseq             |
| In own laboratory    | Commercial kits         | Nextera XT Libray Prep kit (96 samples)***  | NextSeq             |
| In own laboratory    | Commercial kits         | Illumina Nextera XT library Prep Kit  | MiSeq               |

\* 5ng input DNA (as opposed to 1ng)

Altered PCR protocol to favour longer fragment sizes

Adjustment of extension temperature (and final extension) from 72° to 65°C

'Manual' normalisation using library concentration and fragment size as opposed to bead-based normalisation.

\*\* Shearing carried out for 15 minutes at 25°C instead of 20 minutes because 400bp sequencing protocol was used

\*\*\* Half volume for all reagents.

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing

**Table 7. Results of raw reads submitted by participants evaluated by EQA provider QC pipeline summarised by laboratory**

| Parameters                             | Ranges*       | Laboratory ID |           |          |           |          |         |          |           |         |         |         |           |
|--|---------------|---------------|-----------|----------|-----------|----------|---------|----------|-----------|---------|---------|---------|-----------|
|  |               | 19            | 35        | 56       | 70        | 105      | 108     | 129      | 135       | 141     | 142     | 144     | 146       |
| No. of genera detected                 | {1}           | 1             | 1         | 1        | 1         | 1        | 1       | 1        | 1         | 1       | 1       | 1       | 1         |
| Detected species                       | {Lm}          | Lm            | Lm        | Lm       | Lm        | Lm-N     | Lm      | Lm       | Lm        | Lm      | Lm      | Lm      | Lm        |
| Unclassified reads (%)                 |               | 1.5-2.5       | 0.6-2.5   | 0.6-2.2  | 1.5-2.9   | 0.7-50.8 | 1.1-1.8 | 0.6-1.7  | 0.5-1.0   | 0.9-1.8 | 0.8-1.4 | 0.2-1.3 | 0.4-2.0   |
| Length at 25 x min. coverage (Mbp)     | {>2.8 ^ <3.1} | 2.9-3.0       | 2.9-3.0   | 1.8-2.7  | 2.9-3.0   | 0.1-3.0  | 2.9-3.0 | 2.9-2.9  | 2.9-3.0   | 1.0-3.0 | 2.9-3.0 | 2.9-3.0 | 2.9-3.0   |
| Length [0-25] x min. coverage (Mbp)    | {<0.25}       | 0             | 0         | 0        | 0         | 0-0.9    | 0       | 0-0.1    | 0         | 0.0-1.8 | 0       | 0       | 0         |
| No. of contigs at 25 x min. coverage   | {>0}          | 14-21         | 12-25     | 876-1056 | 17-45     | 14-193   | 57-146  | 15-47    | 17-24     | 19-85   | 13-17   | 11-17   | 17-25     |
| No. of contigs [0-25] x min. coverage# | {<1000}       | 0             | 0         | 0        | 0-4       | 0-517    | 0-5     | 0-24     | 0         | 0-165   | 0-2     | 0       | 0-1       |
| Average coverage                       | {>50}         | 160-224       | 40-175    | 61-104   | 51-100    | 8-94     | 30-70   | 50-244   | 153-221   | 24-126  | 40-58   | 75-128  | 140-200   |
| No. of reads (x 1000)                  |               | 1741-2457     | 250-1120  | 707-1278 | 528-1035  | 345-622  | 285-689 | 530-2704 | 1898-2835 | 158-883 | 261-385 | 525-881 | 2148-3169 |
| No. of trimmed reads (x1000)           |               | 1721-2428     | 248-1110  | 691-1235 | 524-1028  | 342-609  | 521-617 | 523-2677 | 1878-2800 | 150-865 | 295-380 | 534-870 | 2148-3169 |
| Maximum read length                    |               | 151           | 301       | 285-365  | 151       | 301      | 241-319 | 151      | 126       | 301     | 251     | 251     | 101       |
| Mean read length                       |               | 140-142       | 215-251   | 217-229  | 143-146   | 204-241  | 186-200 | 139-145  | 123-124   | 218-235 | 245-234 | 210-227 | 97-100    |
| Read insert size                       |               | 267.9-305     | 333-394   | NA       | 288-391   | 199-363  | NA      | 244-450  | 326-351   | 279-358 | 361-399 | 280-327 | 204-360   |
| Insert size StdDev                     |               | 100-106       | 158-199   | NA       | 100-149   | 67-158   | NA      | 108-196  | 175-188   | 102-130 | 157-174 | 93-125  | 85-169    |
| N50 (kbp)                              |               | 238-551       | 274.4-558 | 1.4-3.4  | 162.3-318 | 1.3-407  | 34.0-87 | 125-551  | 295-482   | 22-263  | 262-556 | 353-558 | 286-510   |
| N75 (kbp)                              |               | 143.3-257.3   | 139-263   | 0.9-1.9  | 78-238    | 0.8-262  | 23-45   | 61-258   | 142-258   | 11-236  | 198-262 | 183-263 | 144-262   |

\* Indicative QC range

Lm: *L. monocytogenes*

N: *Neisseria*

# Number of contigs with coverage < 25 (Figure 10B)

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



**Table 4. Results of SNP-based cluster analysis**

| Lab ID   | SNP-based       |  |                          |                          |           |                         |                          |
|----------|-----------------|--|--------------------------|--------------------------|-----------|-------------------------|--------------------------|
|          | Approach        | Reference  | Read mapper              | Variant caller           | Assembler | Distance within cluster | Distance outside cluster |
| Provider | Reference-based | ST6 (REF4)   | BWA                      | GATK                     |           | 0-3                     | 38-71                    |
| 19*      | Reference-based | ST6 ID 2362  | BWA                      | GATK                     |           | 0-4                     | 43-81                    |
| 56       | Assembly-based  |  |                          | ksnp3                    | SPAdes    | 0-57#                   | 561-591 (6109)           |
| 105      | Reference-based | ST6 J1817  | Bowtie2                  | VARSCAN 2                |           | 0-2#                    | 22-42 (1049)             |
| 108      | Reference-based | In-house strain resp ST                                    | CLC assembly cell v4.4.2 | CLC assembly cell v4.4.2 |           | 0-2                     | 37-72                    |
| 142*     | Reference-based | <i>Listeria</i> EGDe (cc9)                                 | CLC Bio                  | CLC Bio                  |           | 0-1219                  | 1223-2814 (8138)         |
| 146      | Reference-based | ST6 ref. CP006046<br>ST1 ref. F2365<br>ST213/ST382 no ref. | BWA                      | In-house                 |           | 0-358                   |                          |

\* Additional analysis

# Only three isolates included due to data quality not meeting laboratory's own QC thresholds

✕ Reported distance to ST6 (non-ST6) isolates (Annex 9).

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



**Table 5. Results of allele-based cluster analysis**


| Lab ID       | Allele based analysis |                             |                           |                               |                           |                            |
|--------------|-----------------------|-----------------------------|---------------------------|-------------------------------|---------------------------|----------------------------|
|              | Approach              | Allelic calling method      | Assembler                 | Scheme                        | Difference within cluster | Difference outside cluster |
| EQA provider | BioNumerics           | Assembly- and mapping-based | SPAdes                    | Applied Math (cgMLST/Pasteur) | 0-3                       | 24-1112                    |
| 19           | BioNumerics           | Assembly- and mapping-based | SPAdes                    | Applied Math (cgMLST/Pasteur) | 0-3                       | 25-1120                    |
| 35           | SeqPhere              | Assembly-based only         | Velvet                    | Ruppitsch (cgMLST)            | 0-2                       | 16-1065                    |
| 70           | SeqPhere              | Assembly-based only         | Velvet                    | Ruppitsch (cgMLST)            | 0-2                       | 16-1062                    |
| 105*         | SeqPhere              | Assembly-based only         | SPAdes v 3.80             | Ruppitsch (cgMLST)            | 0-1 <sup>#</sup>          | 23-812                     |
| 129          | SeqPhere              | Assembly-based only         | Velvet                    | In-house (cgMLST)             | 0-4                       | 15-862                     |
| 135          | SeqPhere              | Assembly-based only         | CLC Genomics Workbench 10 | Ruppitsch (cgMLST)            | 0-2                       | 16-2042                    |
| 141          | SeqPhere              | Assembly-based only         | SPAdes 3.9.0              | Ruppitsch (cgMLST)            | 0-7                       | 19-1060                    |
| 142          | Inhouse               | Assembly-based only         | SPAdes                    | Pasteur (cgMLST)              | 0                         | 13-1120                    |
| 144          | SeqPhere              | Assembly-based only         | Velvet                    | Ruppitsch (cgMLST)            | 0-2                       | 16-1065                    |

\* Additional analysis

<sup>#</sup> Only three isolates included due to data quality not meeting laboratory's own QC thresholds (Annex 9).

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing

**Figure 7. Reported SNP distances or allelic differences for each test isolate to selected cluster representative isolate**



\* Additional analysis

SNP: Single nucleotide polymorphism


Selected cluster representative marked as REF in dark green: Reported cluster of closely related isolates

Light green: Not reported as part of cluster.



# Fifth external quality assessment scheme for *Listeria monocytogenes* typing

## Annex 4. EQA provider cluster analysis based on WGS-derived data



Minimum spanning tree of core genome multi locus sequence typing (cgMLST, [6]) profiles of *L. monocytogenes* EQA-5 isolates.  
Logarithmic scaling in BioNumerics.  
Dark grey: Cluster isolates  
Light grey: Outside cluster isolates

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



## Annex 5. Reported cluster of closely related isolates based on PFGE-derived data

| Lab ID   | Reported cluster                                       | Corresponding REF isolates  | Correct |
|----------|--|---|---------|
| Provider | REF1, REF4, REF9, REF11 (4 and 9 technical duplicates) |   |         |
| 19       | 3362# 2539 2691 2719                                   | REF4, REF1, REF9, REF11   | Yes     |
| 100      | 2080 2295 2405 2499                                    | REF4, REF9, REF11, REF1   | Yes     |
| 105      | 2073 2709 2805 2978                                    | REF4, REF9, REF1, REF11   | Yes     |
| 138      | 2141 2349 2778 2947                                    | REF9, REF1, REF4, REF11   | Yes     |
| 141      | 2022 2050 2092 2872                                    | REF1, REF4, REF11, REF9   | Yes     |
| 142      | 2385 2529 2794 2837                                    | REF9, REF4, REF11, REF1   | Yes     |
| 145      | 2027 2235 2287 2444 2514 2592 2680 2699 2904 2961 2967 | REF5, REF9, REF1, REF3, REF11, REF4<br>REF7, REF2, REF8, REF6, REF 10 | No      |

# Writing error 2362



# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



## Annex 8. Reported cluster of closely related isolates based on WGS-derived data

| Lab ID   | Reported cluster                                       | Corresponding to REF isolates                   | Correct |
|----------|--|---|---------|
| Provider | REF1, REF4, REF9, REF11 (4 and 9 technical duplicates) |   |         |
| 19       | #3562 3539 2691 2719                                   | REF4, REF1 REF9, REF11                          | Yes     |
| 35       | 2251 2737 2783 2993                                    | REF11, REF9, REF1, REF4                         | Yes     |
| 56       | 2341 2165 2612   | REF9, REF1, REF11                               | Yes     |
| 70       | 2104 2216 2567 2767                                    | REF4, REF1, REF11, REF9                         | Yes     |
| 105      | 2073 2805 2978   | REF4, REF1, REF11                               | Yes     |
| 108      | 2098 2788 2582 2422                                    | REF1, REF11, REF9, REF4                         | Yes     |
| 129      | 2079 2640 2912 2950                                    | REF1, REF9, REF11, REF4                         | Yes     |
| 135      | 2161 2423 2673 2897                                    | REF1, REF4, REF11, REF9                         | Yes     |
| 141      | 2022 2050 2092 2872                                    | REF1, REF4, REF11, REF9                         | Yes     |
| 142      | 2385 2529 2794 2837                                    | REF9, REF4, REF11, REF1                         | Yes     |
| 144      | 2143 2626 2727 2822                                    | REF4, REF11, REF1, REF9                         | Yes     |
| 146      | 2068 2197 2377 2488 ##2353 2575 2655 2726              | REF5, REF8, REF3, REF1, REF6, REF4, REF9, REF11 | No      |

#Writing error 2362

##Writing error 2553

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing



## Annex 9. Reported SNP distance and allelic differences

### SNP distances

| Isolate no.        | ST  | Provider | Laboratory ID |      |      |      |      |      |
|--------------------|-----|----------|---------------|------|------|------|------|------|
|                    |     |          | 19*           | 56   | 105  | 108  | 142* | 146  |
| REF1 <sup>†</sup>  | 6   | 3        | 4             | 0*   | 1    | 2    | 1030 | 306  |
| REF2               | 1   | 9999     | 9999          | 9999 | 812  | 9999 | 7502 | 9999 |
| REF3               | 6   | 41       | 47            | 579  | 23   | 45   | 2814 | 0*   |
| REF4 <sup>‡</sup>  | 6   | 0*       | 0*            | 9999 | 1    | 2    | 1219 | 309  |
| REF5               | 6   | 40       | 43            | 561  | 24   | 37   | 2056 | 259  |
| REF6               | 6   | 72       | 81            | 591  | 42   | 72   | 2732 | 358  |
| REF7               | 213 | 9999     | 9999          | 9999 | 734  | 9999 | 8050 | 9999 |
| REF8               | 6   | 39       | 43            | 9999 | 22   | 40   | 1223 | 48   |
| REF9 <sup>‡</sup>  | 6   | 0        | 0             | 57   | 9999 | 2    | 0*   | 296  |
| REF10              | 382 | 9999     | 9999          | 6109 | 745  | 9999 | 8138 | 9999 |
| REF11 <sup>‡</sup> | 6   | 1        | 2             | 48   | 0*   | 0*   | 1114 | 304  |

# Fifth external quality assessment scheme for *Listeria monocytogenes* typing

## Allelic distances

| Isolates no.       | ST  | Provider | Laboratory ID |      |      |      |     |      |      |      |      |
|--------------------|-----|----------|---------------|------|------|------|-----|------|------|------|------|
|                    |     |          | 19            | 35   | 70   | 105* | 129 | 135  | 141  | 142  | 144  |
| REF1 <sup>‡</sup>  | 6   | 3        | 3             | 1    | 1    | 1    | 4   | 1    | 7    | 0    | 1    |
| REF2               | 1   | 1118     | 1120          | 1065 | 1062 | 812  | 862 | 2042 | 1060 | 1120 | 1065 |
| REF3               | 6   | 25       | 25            | 16   | 16   | 23   | 15  | 16   | 19   | 16   | 16   |
| REF4 <sup>‡#</sup> | 6   | 0*       | 0*            | 2    | 2    | 1    | 0   | 2    | 0*   | 0    | 2    |
| REF5               | 6   | 26       | 26            | 18   | 18   | 24   | 17  | 17   | 20   | 13   | 18   |
| REF6               | 6   | 44       | 43            | 35   | 35   | 42   | 32  | 50   | 38   | 41   | 36   |
| REF7               | 213 | 1073     | 1070          | 1028 | 1026 | 734  | 842 | 1031 | 1024 | 1074 | 1028 |
| REF8               | 6   | 28       | 28            | 24   | 24   | 22   | 20  | 39   | 26   | 22   | 24   |
| REF9 <sup>‡#</sup> | 6   | 0        | 0             | 2    | 2    | 9999 | 0*  | 2    | 0    | 0*   | 2    |
| REF10              | 382 | 1060     | 1060          | 1027 | 1021 | 745  | 839 | 1592 | 1025 | 1063 | 1027 |
| REF11 <sup>‡</sup> | 6   | 3        | 2             | 0*   | 0*   | 0*   | 2   | 0*   | 2    | 0    | 0*   |

\* Additional analysis

‡ Closely related isolates

# Technical duplicate isolate

\* Isolate used as cluster representative by participant

9999: Isolates not included in analysis by participant

ST: Sequence type

# Criteria for wg/cgMLST and SNP typing schemes

Schürch et al., CMI, 2018

Examples of relatedness criteria for wg/cgMLST and SNP typing schemes of representative clinically relevant bacteria

| Organism                                     | Relatedness threshold <sup>a</sup>                |      | References  |
|--|---|------|---|
|  | wg/cgMLST (allele)                                | SNPs |   |
| <i>Acinetobacter baumannii</i>               | ≤8  | ≤3   | [25,26]   |
| <i>Brucella</i> spp.                         | Epidemiologic validation in progress <sup>b</sup> |      | <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Campylobacter coli</i> , <i>C. jejuni</i> | ≤14   | ≤15  | [27,28]   |
| <i>Cronobacter</i> spp.                      | Epidemiologic validation in progress <sup>b</sup> |      | <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Clostridium difficile</i>                 | Epidemiologic validation in progress <sup>b</sup> | ≤4   | [29], <a href="http://www.cgmlst.org/ncs">http://www.cgmlst.org/ncs</a> , <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Enterococcus faecium</i>                  | ≤20   | ≤16  | [30]  |
| <i>Enterococcus raffinosus</i>               | Epidemiologic validation in progress <sup>b</sup> |      | <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Escherichia coli</i>                      | ≤10   | ≤10  | [31,32], <a href="https://enterobase.warwick.ac.uk/">https://enterobase.warwick.ac.uk/</a>  |
| <i>Francisella tularensis</i>                | ≤1  | ≤2   | [33,34]   |
| <i>Klebsiella oxytoca</i>                    | Epidemiologic validation in progress <sup>b</sup> |      | <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Klebsiella pneumonia</i>                  | ≤10   | ≤18  | [35,36]   |
| <i>Legionella pneumophila</i>                | ≤4  | ≤15  | [37]  |
| <i>Listeria monocytogenes</i>                | ≤10   | ≤3   | [38,39]   |
| <i>Mycobacterium abscessus</i>               |   | ≤30  | [40]  |
| <i>Mycobacterium tuberculosis</i>            | ≤12   | ≤12  | [41]  |
| <i>Neisseria gonorrhoeae</i>                 | Epidemiologic validation in progress <sup>b</sup> | ≤14  | [42], <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a>   |
| <i>Neisseria meningitidis</i>                | Epidemiologic validation in progress <sup>b</sup> |      | <a href="http://www.cgmlst.org/ncs">http://www.cgmlst.org/ncs</a>   |
| <i>Pseudomonas aeruginosa</i>                | ≤14   | ≤37  | [31,43]   |
| <i>Salmonella dublin</i>                     | Epidemiologic validation in progress <sup>b</sup> | ≤13  | [44], <a href="https://enterobase.warwick.ac.uk/">https://enterobase.warwick.ac.uk/</a>   |
| <i>Salmonella enterica</i>                   | Epidemiologic validation in progress <sup>b</sup> | ≤4   | [45], <a href="http://www.cgmlst.org/ncs">http://www.cgmlst.org/ncs</a> , <a href="http://www.applied-maths.com/applications/wgmlst">http://www.applied-maths.com/applications/wgmlst</a> , <a href="https://enterobase.warwick.ac.uk/">https://enterobase.warwick.ac.uk/</a> |
| <i>Salmonella typhimurium</i>                | Epidemiologic validation in progress <sup>b</sup> | ≤2   | [46], <a href="https://enterobase.warwick.ac.uk/">https://enterobase.warwick.ac.uk/</a>   |
| <i>Staphylococcus aureus</i>                 | ≤24   | ≤15  | [47,48]   |
| <i>Streptococcus suis</i>                    |   | ≤21  | [49]  |
| <i>Vibrio parahaemolyticus</i>               | ≤10   |      | [50]  |
| <i>Yersinia</i> spp.                         | 0   |      | [51]  |

# Thanks for your attention!

---